

ARG82910 Human Toxoplasma IgG antibody ELISA Kit

Package: 96 wells

Store at: 4°C

Summary

Product Description	ARG82910 Human Toxoplasma IgG antibody ELISA Kit has been designed for the qualitative determination of specific IgG antibodies against Toxoplasma in serum and plasma.
Tested Reactivity	Hu
Tested Application	ELISA
Target Name	Toxoplasma IgG antibody
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 453 nm
Sensitivity	Diagnostic Sensitivity: 100 %
Sample Type	Serum and plasma (Citrate, heparin).
Standard Range	Cut - off
Sample Volume	5-10 µl

Application Instructions

Assay Time	1 h, 30, 15 min
------------	-----------------

Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	for laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Highlight	Related products: Toxoplasma ELISA Kits; New ELISA data calculation tool: Simplify the ELISA analysis by GainData
Research Area	Microbiology and Infectious Disease kit