


## ARG82986 Human IL10 ELISA Kit

Package: 96 wells  
Store at: 4°C or -20°C  
(depending on the component)

### Summary

Product Description	ARG82986 Human IL10 ELISA Kit is an Enzyme Immunoassay kit for the quantification of Human IL10 in Serum, Plasma, Cell culture supernatants.
Tested Reactivity	Hu
Tested Application	ELISA
Target Name	IL10
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 450 nm
Sensitivity	4 pg/ml
Sample Type	Serum, Plasma, Cell culture supernatants
Standard Range	7.8 - 500 pg/ml
Sample Volume	50 µl
Precision	

### Images



ARG82986 Human IL10 ELISA Kit standard curve image

ARG82986 Human IL10 ELISA Kit results of a typical standard run with optical density reading at 450 nm.