

## ARG83000 Hyaluronidase Inhibitor Screening Assay Kit (Colorimetric)

Package: 100 tests  
Store at: -20°C

### Summary

Product Description	ARG83000 Hyaluronidase Inhibitor Screening Assay Kit is a detection kit for the quantification of Hyaluronidase Inhibitor Screening in Test Compounds, Hyaluronidase Inhibitors samples.
Tested Reactivity	Other
Tested Application	FuncSt
Target Name	Hyaluronidase Inhibitor Screening
Conjugation Note	Read at 600 nm (colorimetric)
Sample Type	Test Compounds, Hyaluronidase Inhibitors samples.
Sample Volume	20 µl

### Application Instructions


Application Note	Please note that this kit does not include a microplate.
Assay Time	45 mins

### Properties

Form	Liquid
Storage instruction	Store the kit at -20°C. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

### Bioinformation

Highlight	Related news: <a href="#">TCA intermediate fumarate promotes mitobiogenesis</a>
-----------	--



ARG83000 Hyaluronidase Inhibitor Screening Assay Kit (Colorimetric) enzyme activity assay image

ARG83000 Hyaluronidase Inhibitor Screening Assay Kit (Colorimetric) results of a typical example run with optical density reading at 660 nm.