

ARG83071
Feline SARS-CoV-2 antibody ELISA Kit

Package: 96 wells

Store at: 4°C

Summary

Product Description	ARG83071 Feline SARS-CoV-2 antibody ELISA Kit is an Enzyme Immunoassay kit for qualitative of SARS-CoV-2 antibody in cat serum.
Tested Reactivity	Cat
Tested Application	ELISA
Target Name	Feline SARS-CoV-2
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 450 nm.
Detection Range	Cut - off
Sample Type	Serum
Sample Volume	100 µl
Precision	Intra-Assay CV: less than 5% Inter-Assay CV: less than 8%
Alternate Names	SARS-CoV

Application Instructions

Assay Time	~2 hour
------------	---------

Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Background	Severe acute respiratory syndrome coronavirus 2 (SARS-CoV-2) is a strain of coronavirus that causes COVID-19 (coronavirus disease 2019), the respiratory illness responsible for the ongoing COVID-19 pandemic. The virus previously had a provisional name, 2019 novel coronavirus (2019-nCoV), and has also been called the human coronavirus 2019 (HCoV-19 or hCoV-19). First identified in the city of Wuhan, Hubei, China, the World Health Organization declared the outbreak a public health emergency of international concern on January 30, 2020, and a pandemic on March 11, 2020. SARS-CoV-2 is a positive-sense single-stranded RNA virus that is contagious in humans.
Highlight	Related products: SARS-CoV antibodies SARS-CoV ELISA Kits SARS-CoV recombinant proteins;

Related news:

[HMGB1, a biomarker and therapeutic target in COVID-19](#)

[Exploring Antiviral Immune Response](#)

New ELISA data calculation tool:

[Simplify the ELISA analysis by GainData](#)