

## ARG83078 FSME (TBEV) ELISA Kit

Package: 96 wells  
Store at: 4°C

### Summary

Product Description	ARG83078 FSME (TBEV) ELISA Kit is an Enzyme Immunoassay kit for the qualitative of FSME (TBEV) in animal serum.
Tested Reactivity	Other
Tested Application	ELISA
Target Name	FSME (TBEV)
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 450 nm.
Detection Range	Cut - off
Sample Type	Serum
Sample Volume	100 µl
Precision	Intra-Assay CV: less than 4 % Inter-Assay CV: less than 10 %
Alternate Names	Tick borne encephalitis virus (TBEV)

### Application Instructions

Assay Time	~2 hour
------------	---------

### Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

### Bioinformation

Highlight	Related products: <a href="#">FSME ELISA Kits;</a> New ELISA data calculation tool: <a href="#">Simplify the ELISA analysis by GainData</a>
-----------	--