

ARG83092 E.coli Residual Total RNA Detection Kit

Package: 100 tests
Store at: -20°C

Summary

Product Description	E.coli Residual Total RNA Detection Kit is designed for the quantitative detection of residual E.coli total RNA in various biological products to improve control quality of nucleic acid. This kit adopts the principle of the RT-PCR fluorescent probe, combining reverse transcription PCR technology and fluorescent probe method, to realize one-step quantitative detection.
Tested Reactivity	E. coli
Tested Application	FuncSt
Target Name	E.coli Residual Total RNA
Sensitivity	5×10^{-4} pg/ μ L
Standard Range	2×10^{-3} - 2×10^2 pg/ μ L

Properties

Storage instruction	Store the kit at -20°C at all times. All reagents must be kept on ice during the entire experiment.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Highlight	Related products: Quality control of transgene DNA; Related news: Quality control for engineered T cell therapy
-----------	--