

ARG83106 RNase Inhibitor ELISA Kit

Package: 96 tests
Store at: 4°C

Summary

Product Description	RNase Inhibitor ELISA Kit is designed for the quantitative detection of residual RNase Inhibitor content added in RNA pharmaceuticals processes by using a double-antibody sandwich method.
Tested Application	ELISA
Target Name	RNase Inhibitor
Conjugation	HRP
Sensitivity	5 ng/mL
Standard Range	20 - 640 ng/mL


Application Instructions

Assay Time	~3 hours
------------	----------

Properties

Storage instruction	Store the kit at 4°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Images



ARG83106 RNase Inhibitor ELISA Kit standard curve image

ARG83106 RNase Inhibitor ELISA Kit results of standard run with optical density reading at 450 nm