

ARG83507 Feline IgE ELISA Kit

Package: 96 wells
Store at: 4°C

Summary

Product Description	ARG83507 Feline IgE ELISA Kit is an Enzyme Immunoassay kit for the quantification of Feline IgE in Plasma and Serum
Tested Reactivity	Cat
Tested Application	ELISA
Target Name	IgE
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 450 nm.
Sensitivity	6.25 ng/ml
Sample Type	Plasma and Serum
Standard Range	12.5 - 400 ng/ml
Sample Volume	100 µl
Precision	Intra-Assay CV: less than 10% Inter-Assay CV: less than 10%
Alternate Names	IgE; Immunoglobulin E

Application Instructions

Assay Time	2.5 hours
------------	-----------

Properties

Form	96 well
Storage instruction	Store the kit at 4°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Gene Symbol	IgE
Gene Full Name	Immunoglobulin E