

## ARG83564 Histamine Dehydrogenase Activity Assay Kit

Package: 100 tests

Store at: 4°C

### Summary

Product Description	ARG83564 Histamine Dehydrogenase Activity Assay Kit can be used to measure Histamine Dehydrogenase in Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
Tested Reactivity	Other
Tested Application	FuncSt
Target Name	Histamine Dehydrogenase
Conjugation Note	Read at 450 nm
Sensitivity	0.005 µmol/ml
Sample Type	Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
Standard Range	0.01 - 1 µmol/ml
Sample Volume	10 µl
Alternate Names	EC 1.4.99

### Application Instructions

Assay Time	5 min
------------	-------

### Properties

**Storage instruction** Store the kit at 4°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.

**Note** For laboratory research only, not for drug, diagnostic or other use.