

ARG83566 beta HB Assay Kit

Package: 100 tests
Store at: 4°C

Summary

| | |
|---------------------|---|
| Product Description | ARG83566 beta HB Assay Kit can be used to measure beta HB in Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids |
| Tested Reactivity | Other |
| Tested Application | FuncSt |
| Target Name | beta HB |
| Conjugation Note | Read at 450 nm |
| Sensitivity | 0.025 µmol/ml |
| Sample Type | Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids |
| Standard Range | 0.05 - 5 µmol/ml |
| Sample Volume | 10 µl |

Application Instructions

| | |
|------------|-------|
| Assay Time | 5 min |
|------------|-------|

Properties

| | |
|---------------------|--|
| Storage instruction | Store the kit at 4°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components. |
| Note | For laboratory research only, not for drug, diagnostic or other use. |