

ARG83607 2-Phosphoglycerate Assay Kit

Package: 100 assay
Store at: 4°C, -20°C

Summary

Product Description	ARG83607 2-Phosphoglycerate Assay Kit is a detection kit for the quantification of 2-Phosphoglycerate.
Tested Application	FuncSt
Target Name	2-Phosphoglycerate
Conjugation Note	Read at 540 nm.
Sensitivity	1.57 μ M
Sample Type	Serum, plasma, urine and cell culture supernatants.
Standard Range	3.13 - 200 μ M
Sample Volume	50 μ L
Alternate Names	3-hydroxy-2-phosphonoxypropanoic acid; 2-Phosphoglyceric acid; 2PG

Properties

Form	96 well
Storage instruction	Store the kit at 4°C, -20°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.